

## Phospho-BACE1 (Thr252) Antibody

Cat.#: AF3650	Concn.: 1mg/ml	Mol.Wt.:
Size: 100ul,200ul	Source: Rabbit	Clonality: Polyclonal

Application:	ELISA(peptide) 1:20000-1:40000 *The optimal dilutions should be determined by the end user.
Reactivity:	Human,Mouse,Rat
Purification:	The antibody is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.
Specificity:	Phospho-BACE1 (Thr252) Antibody detects endogenous levels of BACE1 only when phosphorylated at Thr252.
Immunogen:	A synthesized peptide derived from human BACE1 around the phosphorylation site of Thr252.
Uniprot:	P56817
Storage Condition and Buffer:	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.Store at -20 °C.Stable for 12 months from date of receipt.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.