

## **Goat Anti-Mouse IgG (H+L) Fluor594-conjugated**

Cat.#: S0005  
Size: 50ul,100ul,200ul

Concn.: 1mg/ml  
Source: Goat

Mol.Wt.:  
Clonality: Polyclonal

Application: WB 1:100-1:500, IHC 1:50-1:200, IF/ICC 1:100-1:500

Reactivity: Mouse

Purification: affinity purification.

Immunogen: The antiserum was developed in goat using the MOUSE IgG as the immunogen.

Storage Condition and Buffer: 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. Store at -20°C□avoid repeated thawing and freezing.

**IMPORTANT:** For western blot, incubate membrane with diluted secondary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 37°C, 2H.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.