

Goat Anti-Mouse IgG (H+L) Fluor488-conjugated

Cat.#: S0017
Size: 100ul,200ul

Concn.: 1mg/ml
Source: Goat

Mol.Wt.:
Clonality: Polyclonal


Application: WB 1:100 - 1:500, IF/ICC 1:100 - 1:500

Reactivity: Mouse

Purification: affinity purification.

Immunogen: The antiserum was developed in goat using the MOUSE IgG as the immunogen.

Storage Condition and Buffer: 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. Store at -20°C avoid repeated thawing and freezing.



immunofluorescence analysis of HeLa cells using tubulin mouse mAb.
Green: Goat Anti-Mouse IgG(H+L) Alexa Fluor488 -conjugated

IMPORTANT: For western blot, incubate membrane with diluted secondary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 37°C, 2H.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.