

Strep tag II antibody

Cat.#: T0056
 Size:

Concn.: 1mg/ml
 Source: Mouse

Mol.Wt.:
 Clonality: Monoclonal

Application: WB(1:1000-1:10000)

Reactivity: All

Purification: Protein G purification.


Storage Condition and Buffer: Mouse IgG1 in 10 mM PBS (pH 7.4), 0.1mg/ml BSA , 50% glycerol and 1% ProClin™ 300.Store at -20 °C.Stable for 12 months from date of receipt.

Strep tag II fusion protein



(Strep tag II)1:5000

Western blot analysis of strep tag fusion protein, using Strep tag II mouse monoclonal antibody.



Western blot- Strep tag II mouse monoclonal antibody

Western blot analysis of strep tag fusion protein, using Strep tag II mouse monoclonal antibody.

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.