


HPLC REPORT

Sample: G-13 Analyzed date: 2014-10-11
Analyst: HXH
Lot. No.: P140922-WF608838
Column: SinoChrom ODS-BP, 4.6*250mm, 5µm
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
Gradient:
0.0min 15% 85%
25.0min 40% 60%
25.1min 100% 0%
30.0min Stop
Volume: 20µl
Wavelength: 220nm
Flow rate: 1.0ml/min



Rank	Time	Area	Conc.
1	8.255	68431.211	2.1955
2	9.075	72472.445	2.3252
3	9.562	93735.445	3.0073
4	9.828	2812949.000	90.2486
5	9.828	69300.523	2.2234
Total		3116888.625	100