


# MASS SPECTROMETRY REPORT



## Sample Description

Analyzed date: 2017-05-14  
 Analyst: YU  
 Sample: C-12  
 M.W.: 1290.38  
 Lot. No.: P170509-TL654058

## Instrument

Agilent-6125B  
 Probe: ESI  
 Nebulizer Gas Flow: 1.5L/min  
 CDL: -20.0v  
 CDL Temp.: 250 °C  
 Block Temp.: 200 °C

Probe Bias: + 4.5kv  
 Detector: 1.5kv  
 T. Flow: 0.2ml/min  
 B. Conc.: 50% H2O/50% ACN