

Phospho-EPHA6 (Tyr978) Ab

[Images\(1\)](#)

Cat.#: AF8427	Concn.: ~1mg/ml	Mol.Wt.: 120~150kDa
Size:	Source: Rabbit	Clonality: Polyclonal

Application: WB 1:1000-3000

*The optimal dilutions should be determined by the end user.

Reactivity: Human, Mouse, Rat

Storage: Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification: The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.

Immunogen: A synthesized peptide derived from human EPHA6 around the phosphorylation site of Tyr978.

Uniprot: Q9UF33

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.