


## CD52 Ab

[Images\(4\)](#)


|               |                 |                       |
|---------------|-----------------|-----------------------|
| Cat.#: DF8797 | Concn.: ~1mg/ml | Mol.Wt.: 7 kDa, 20kDa |
| Size:         | Source: Rabbit  | Clonality: Polyclonal |

|               |                                                                                                                                                                |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application:  | WB 1:1000-3000, IHC 1:50-1:200<br>*The optimal dilutions should be determined by the end user.                                                                 |
| Reactivity:   | Human, Monkey                                                                                                                                                  |
| Storage:      | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |
| Purification: | The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).                                      |
| Immunogen:    | A synthesized peptide derived from human CD52, corresponding to a region within the internal amino acids.                                                      |
| Uniprot:      | P31358                                                                                                                                                         |



Western blot analysis of extracts from VERO cells, using CD52 Ab. The lane on the left was treated with blocking peptide.

Observed bands: 20 kDa.



DF8797 at 1/100 staining Human gastric cancer by IHC-P. The sample was formaldehyde fixed and a heat mediated antigen retrieval step in citrate buffer was performed. The sample was then blocked and incubated with the primary Ab at 4°C overnight. An HRP conjugated anti-Rabbit Ab was used as the secondary Ab.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.