



## FHL3 Ab

[Images\(3\)](#)

|               |                                                                                                                                                               |                       |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|
| Cat.#: DF8826 | Concn.: ~1mg/ml                                                                                                                                               | Mol.Wt.: 31 kDa       |
| Size:         | Source: Rabbit                                                                                                                                                | Clonality: Polyclonal |
| Application:  | IHC 1:50-1:200, WB 1:1000-3000<br>*The optimal dilutions should be determined by the end user.                                                                |                       |
| Reactivity:   | Human, Mouse, Rat                                                                                                                                             |                       |
| Storage:      | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |                       |
| Purification: | The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).                                     |                       |
| Immunogen:    | A synthesized peptide derived from human FHL3, corresponding to a region within the internal amino acids.                                                     |                       |
| Uniprot:      | Q13643                                                                                                                                                        |                       |



Western blot analysis of extracts from Rat heart, using FHL3 Ab. The lane on the left was treated with blocking peptide.



DF8826 at 1/100 staining rat heart tissue by IHC-P. The sample was formaldehyde fixed and a heat mediated antigen retrieval step in citrate buffer was performed. The sample was then blocked and incubated with the primary Ab at 4°C overnight. An HRP conjugated anti-rabbit Ab was used as the secondary Ab.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.