


## KCNJ10 Ab


[References\(1\)](#) [Images\(6\)](#)

|               |                 |                       |
|---------------|-----------------|-----------------------|
| Cat.#: DF9260 | Concn.: ~1mg/ml | Mol.Wt.: 43 kDa       |
| Size:         | Source: Rabbit  | Clonality: Polyclonal |

|               |                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application:  | WB 1:1000-3000, IHC 1:50-1:200<br>*The optimal dilutions should be determined by the end user.                                                                |
| Reactivity:   | Human, Mouse, Rat                                                                                                                                             |
| Storage:      | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |
| Purification: | The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).                                     |
| Immunogen:    | A synthesized peptide derived from human KCNJ10, corresponding to a region within the internal amino acids.                                                   |
| Uniprot:      | P78508                                                                                                                                                        |



Western blot analysis of extracts from HepG2, using KCNJ10 Ab. Lane 1 was treated with the blocking peptide.



DF9260 at 1/100 staining Rat colorectal tissue by IHC-P. The sample was formaldehyde fixed and a heat mediated antigen retrieval step in citrate buffer was performed. The sample was then blocked and incubated with the primary Ab at 4°C overnight. An HRP conjugated anti-Rabbit Ab was used as the secondary Ab.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.