

Goat Anti-Mouse IgG (H+L) AP

[References\(2\)](#)

Cat.#: S0004
Size:

Concn.: 1mg/ml
Source: Goat

Mol.Wt.:
Clonality: Polyclonal

Application:

WB: 1:5000~1:10000

*The optimal dilutions should be determined by the end user.

Reactivity:

Mouse

Storage:

AP-labeled Goat IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl. Store at -20°C, avoid repeated thawing and freezing.

Purification:

affinity purification.

Immunogen:

The antiserum was developed in goat using the mouse IgG as the immunogen.

Conjugate:

AP

IMPORTANT: For western blot, incubate membrane with diluted secondary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 37°C, 2H.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.