

RFP-Tag Ab

[References\(3\)](#) [Images\(3\)](#)

Cat.#: T0055
Size:

Concn.: ~1mg/ml
Source: Mouse

Mol.Wt.:
Clonality: Monoclonal

Application:

WB(1:1000-1:10000)

*The optimal dilutions should be determined by the end user.

Reactivity:


All

Storage:

Mouse IgG1 in 10 mM PBS (pH 7.4), 0.1mg/ml BSA , 50% glycerol and 1% ProClin™ 300. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:

Protein G purification.



Western blot analysis of RFP Mouse Monoclonal Ab expression in RFP-tag fusion protein.

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.