

## Strep tag II Ab

[Images\(2\)](#)

Cat.#: T0056  
Size:

Concn.: ~1mg/ml  
Source: Mouse

Mol.Wt.:  
Clonality: Monoclonal

Application:

WB(1:1000-1:10000)

\*The optimal dilutions should be determined by the end user.

Reactivity:

All

Storage:

Mouse IgG1 in 10 mM PBS (pH 7.4), 0.1mg/ml BSA , 50% glycerol and 1% ProClin™ 300. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:

Protein G purification.

Strep tag II  
fusion protein

Affinit

(Strep tag II)1:5000

Western blot analysis of strep tag fusion protein, using Strep tag II mouse monoclonal Ab.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.