

EBFP Ab

[Images\(1\)](#)

Cat.#: BF0706	Concn.: ~1mg/ml	Mol.Wt.:
Size:	Source: Mouse	Clonality: Monoclonal

Application: WB 1:2000-1:10000
*The optimal dilutions should be determined by the end user.

Reactivity: All

Storage: Mouse IgG1 in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification: Affinity-chromatography.

Immunogen: Recombinant EBFP expressed in E. coli.



Western blot analysis of recombinant EBFP proteins, using EBFP mouse monoclonal Ab.

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.