

TNNI3 Ab

[Images\(1\)](#)

Cat.#: AF6215
Size:

Concn.: ~1mg/ml
Source: Rabbit

Mol.Wt.: 26kDa
Clonality: Polyclonal

Application:

WB 1:500-1:2000, IHC 1:50-1:200
*The optimal dilutions should be determined by the end user.

Reactivity:

Human, Mouse, Rat

Storage:

Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.

Purification:

The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).

Immunogen:

A synthesized peptide derived from human TNNI3, corresponding to a region within N-terminal amino acids.

Uniprot:

P19429

Description:

Troponin I (TnI) is one of 3 subunits that form the troponin complex of the thin filaments of striated muscle; see also TNNI1 (MIM 191042). The others are troponin T (TnT; see MIM 191041) and troponin C (TnC; see MIM 191040).

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.