

FGFP1 Ab

[References\(1\)](#) [Images\(2\)](#)

| | | |
|---------------|-----------------|-----------------------|
| Cat.#: DF2498 | Concn.: ~1mg/ml | Mol.Wt.: 26-38kD |
| Size: | Source: Rabbit | Clonality: Polyclonal |

| | |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application: | WB 1:500-1:2000, IHC 1:50-1:200 *The optimal dilutions should be determined by the end user. |
| Reactivity: | Human |
| Storage: | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |
| Purification: | The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific). |
| Immunogen: | A synthesized peptide derived from human FGFP1, corresponding to a region within C-terminal amino acids. |
| Uniprot: | Q14512 |

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.