


## HPLC REPORT

Sample: C-9  
Analyst: HXH  
Lot. No.: P140710-WF594431  
Column: SinoChrom ODS-BP, 4.6\*250mm, 5µm  
Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile  
Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water  
Gradient:  
0.0min 10% 90%  
25.0min 35% 65%  
25.1min 100% 0%  
30.0min Stop  
Volume: 20µl  
Wavelength: 220nm  
Flow rate: 1.0ml/min

Analyzed date: 2014-07-13



| Rank  | Time  | Area        | Conc.   |
|-------|-------|-------------|---------|
| 1     | 6.650 | 5254.396    | 0.2504  |
| 2     | 8.768 | 26233.344   | 1.2503  |
| 3     | 9.097 | 70237.266   | 3.3475  |
| 4     | 9.260 | 1862649.375 | 88.7749 |
| 5     | 9.260 | 133797.047  | 6.3768  |
| Total |       | 2098171.427 | 100     |