


# MASS SPECTROMETRY REPORT



Sample Description

Analyzed date: 2014-3-8

Analyst: SHAO

Sample: A-13

M.W.: 1491.80

Lot. No.: P140306-SY566575

Instrument

SHIMADZU LCMS-2020

Probe:

ESI

Probe Bias:

+ 4.5kv

Nebulizer Gas Flow:

1.5L/min

Detector:

1.5kv

CDL:

-20.0v

T. Flow:

0.2ml/min

CDL Temp.:

250°C

B. Conc.:

50% H<sub>2</sub>O/50% ACN

Block Temp.:

200°C