


HPLC REPORT

Sample: D-15 Analyzed date: 2016-05-18
 Analyst: LY
 Lot. No.: P160509-YS654265
 Column: Kromasil-C18, 4.6*250mm, 5µm
 Solvent A: A: 0.1% Trifluoroacetic Acid in 100% Acetonitrile
 Solvent B: B: 0.1% Trifluoroacetic Acid in 100% Water
 Gradient:

	A	B
0.0min	22%	78%
25.0min	47%	53%
25.1min	100%	0%
30.0min	Stop	

Volume: 10µl
 Wavelength: 220nm
 Flow rate: 1.0ml/min



Peak No	Ret Time	Height	Area	Conc.
1	8.878	2752.501	31820.254	0.2955
2	9.397	7537.286	117123.094	1.0878
3	9.747	11298.726	131492.969	1.2212
4	9.932	931033.938	9764343.000	90.6864
5	9.932	98405.992	722375.438	6.7091
Total				100.00