

KT3-Tag Ab

[Images\(1\)](#)

Cat.#: T0011	Concn.: ~1mg/ml	Mol.Wt.:
Size:	Source: Mouse	Clonality: Monoclonal

Application:	WB 1:3000-1:10000 *The optimal dilutions should be determined by the end user.
Reactivity:	All
Storage:	Mouse IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt.
Purification:	Affinity-chromatography.
Immunogen:	A synthetic peptide KPPTPPPEPET coupled to KLH.
Description:	KT3-tag peptide KPPTPPPEPET conjugated to KLH. KT3-Tag antibody can recognize C-terminal, internal, and N-terminal KT3-tagged proteins.



1µg KT3 fusion protein+ Primary Ab dilution at 1?1:5,000 2?1:10,000

IMPORTANT: For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.