

## Phospho-TrkC (Tyr516) Ab

|               |  |                       |
|---------------|--|-----------------------|
| Cat.#: AF3582 | Concn.: ~1mg/ml  | Mol.Wt.:              |
| Size:         | Source: Rabbit   | Clonality: Polyclonal |
| Application:  | ELISA(peptide) 1:20000-1:40000<br>*The optimal dilutions should be determined by the end user.   |                       |
| Reactivity:   | Human,Mouse,Rat  |                       |
| Storage:      | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |                       |
| Purification: | The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns.       |                       |
| Immunogen:    | A synthesized peptide derived from human TrkC around the phosphorylation site of Tyr516.   |                       |
| Uniprot:      | Q16288   |                       |

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.