


## p23 Ab

[Images\(1\)](#)

|               |                 |                       |
|---------------|-----------------|-----------------------|
| Cat.#: AF6794 | Concn.: ~1mg/ml | Mol.Wt.:              |
| Size:         | Source: Rabbit  | Clonality: Polyclonal |

|               |                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application:  | WB 1:500-1:2000<br>*The optimal dilutions should be determined by the end user.                                                                               |
| Reactivity:   | Human, Mouse, Rat                                                                                                                                             |
| Storage:      | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 °C. Stable for 12 months from date of receipt. |
| Purification: | The antiserum was purified by peptide affinity chromatography using SulfoLink™ Coupling Resin (Thermo Fisher Scientific).                                     |
| Immunogen:    | A synthesized peptide derived from human p23.                                                                                                                 |
| Uniprot:      | Q15185                                                                                                                                                        |



Western blot analysis of extracts from MCF7 cells (serum starvation treatment), using p23 Ab. The lane on the left was treated with blocking peptide.

**IMPORTANT:** For western blot, incubate membrane with diluted primary Ab in 5% w/v milk, 1X TBS, 0.1% Tween@20 at 4°C with gentle shaking, overnight.

For Research Use Only. Not for use in diagnostic and therapeutic procedures. Not for resale without express authorization.